



Reference: PTXBC007451

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPBD1A; MBDIN; NTPBP; RPAP4; XAB1; GPN-loop GTPase 1; ATP(GTP)-binding protein; MBD2-interacting protein; RNA polymerase II associated protein 4; RNAPII-associated protein 4; XPA binding protein 1, GTPase		

For research use only.

GENE DESCRIPTION	GPN-loop GTPase 1
GENE ID	11321
ORIGIN SPECIES	Human
ACCESSIONS	BC007451
VECTOR	pDONR223
NCBI SYMBOL	GPN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]