

OSTF1-osteoclast stimulating factor 1 Gene

Reference: PTXBC007459

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OSF; SH3P2; bA235O14.1; osteoclast-stimulating factor 1; osteoclast stimulating factor 1		

For research use only.

GENE DESCRIPTION	osteoclast stimulating factor 1
GENE ID	26578
ORIGIN SPECIES	Human
ACCESSIONS	BC007459
VECTOR	pDONR223
NCBI SYMBOL	OSTF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaagccgcccccaaccagtcgaaccaggtgagggagggaagttaaagtcctcagagccctgtatagctgttgaaaccagaactccagatgaattatactttgaggaaggtgatattatctacattactgacatgagcgataccaattggtgaaaggcacctccaaaggcaggac
tggactaatccaagcaactatgtgctgagcaggcagaatccattgacaatccattgcatgaagcagcaaaaaggcaactgagctggttgagagaggttttgacacagagtggtgttaatggcttagacaaagctggaagcactgccttatactgggttgccacggggccac
aaagatatatggaaatgctattactcaacaaatattgaactgaaccagcagaacaagttgggagatacagcttgcagctgctgcctggaagggttatgcagatcgtccagttgttctggcaaaagggtctagaacagacttaagaacattgagaagaagctggccttcgacatgg
ctaccaatgctgctgcatctctctgaaaaaagaacagggaacagatgcagttcgaacattaagcaatgccgaggactatctcgtgatgaagactcagattaa