

NSUN5C-NOL1/NOP2/Sun domain family, member 5C Gene

Reference: PTXBC007488

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	NOL1/NOP2/Sun domain family, member 5C
GENE ID	260294
ORIGIN SPECIES	Human
ACCESSIONS	BC007488
VECTOR	pDONR223
NCBI SYMBOL	NSUN5C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggagctgctgtgtttcccccagacagatctgcatgaacacccactgtacggggccggacacctcattctgcaggacagggccagctgtctccagccatgctgtggaccccgccaggtcccatgtcatggatgcctgtgccacccaggaataaagaccagtcacttg
gctgctcttctgaagaaccaagggaagatcttgcctttgacctggatgccaggcggctggcatccatggccacgctgctggcctgggctggcgtctctctgtgtgagctggctgaggaggacttctggcggtctcccccttagatccgcgctatcgtgaggtccactatgctgctggat
ccttctgcagtggtcgggtgagatggtatgccgagcagacagctggaggagccggggcaggggacacctaagccgggtgctgtcatgcccctggcagggttccagcagcgagccctgtgccacgcgctcacttccctccctgcagcgggtcgtctactccatgtgctccctgcc
aggaggagaatgaagacatgtacaagatgcgtgcagcagaacccggcgcttcaggctag