



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

UBE2S-ubiquitin-conjugating enzyme E2S Gene

Reference: PTXBC007554

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

E2-EPF; E2EPF; EPF5; ubiquitin-conjugating enzyme E2 S; E2 ubiquitin-conjugating enzyme S; E2-EPF5; ubiquitin carrier protein S; ubiquitin conjugating enzyme E2S; ubiquitin-conjugating enzyme E2-24 kD; ubiquitin-conjugating enzyme E2-24 kDa; ubiquitin-con

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2S
GENE ID	27338
ORIGIN SPECIES	Human
ACCESSIONS	BC007554
VECTOR	pDONR223
NCBI SYMBOL	UBE2S
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr