

ATG4D-ATG4 autophagy related 4 homolog D (*S. cerevisiae*) Gene

Reference: PTXBC007639

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cysteine protease ATG4D; APG4-D; APG4D; AUTL4; APG4 autophagy 4 homolog D; ATG4 autophagy related 4 homolog D; AUT-like 4 cysteine endopeptidase; autophagin-4; autophagy-related cysteine endopeptidase 4; autophagy-related protein 4 homolog D; cysteine pro		

For research use only.

GENE DESCRIPTION	ATG4 autophagy related 4 homolog D (<i>S. cerevisiae</i>)
GENE ID	84971
ORIGIN SPECIES	Human
ACCESSIONS	BC007639
VECTOR	pDONR223
NCBI SYMBOL	ATG4D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgggaaacgcgacactcactgtactcattgggtaccaagatgacttctgctgtacctggaccctcactactgccagccactgtggatgtcagccagccgactccccctggagtcctccactgcacctcgccccgcaagatggcctttgccaagatggaccaagctgtacc
gtgggcttctatgctggagacaggaaggagtttgagacactctgtcagagctgaccagggtctcctcagctcctcagccagagcggtaccccatgttcacacctggccgaggccatgctcaggaccacagcctggacgacctctgtccccagctcgccccgccacactccgct
ccctgcacagggcggtcctcagggccaacgccccagctctgaggactttgttttataa