

RBM14-RNA binding motif protein 14 Gene

Reference: PTXBC007641

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COAA; PSP2; SIP; SYTIP1; TMEM137; RNA-binding protein 14; RRM-containing coactivator activator/modulator; SYT-interacting protein; paraspeckle protein 2; synaptotagmin-interacting protein; transmembrane protein 137; RNA binding motif protein 14		

For research use only.

GENE DESCRIPTION	RNA binding motif protein 14
GENE ID	10432
ORIGIN SPECIES	Human
ACCESSIONS	BC007641
VECTOR	pDONR223
NCBI SYMBOL	RBM14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagaacttactgatgctcttttgggaaaaagggtgctttaagactcgattgggagcctactggggcaaggcgactggggaggtgtacgtgttagtcaagtggaggtcaggtgggttatctctgcctcctcactggctcatgtgaactgctgtgacagaatatcttgcctagatgcc
tcttccctcttgccatcgtgcaaaagcattcgatcttatga