



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ARMCX6-armadillo repeat containing, X-linked 6 Gene

Reference: PTXBC007677

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein ARMCX6; GASP10; armadillo repeat containing, X-linked 6		
For research use only.			

GENE DESCRIPTION	armadillo repeat containing, X-linked 6
GENE ID	54470
ORIGIN SPECIES	Human
ACCESSIONS	BC007677
VECTOR	pDONR223
NCBI SYMBOL	ARMCX6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcggcgtcgggaagtgggtggatgctggcaggactgatgattggggctgtgctgctactgcgtttacaactgaccatagaaagatgacagtggagaagctggaggaggagggggaagagaggtgtggcagatgaccacggagctggatgaggaggagcgtgatattg
gttgattttgaacctatgctcggccctgggactgaggatggggattgggactgaacctggggcccccagttggcactgaggacagccctcagggggaggcaaggcccaaccgagcaccccaataaaacagcggccattcccctatgaacataaaatacttggagtgctcaaaattgta
aaaatggcactgtgtcttgcaactcccaagtgtcttttcacagggaacaactgttttgcgcagccgaagatcgccgcttccatttggcagcagatatacatgaatgctgtgcagcctctcaatgctgataaacacagccccactgcagaccctctgtagaagagcgtttgtgtgcc
cgcagaacttaaatgcgagtattgaagctggagcattgaagatgcatcatgttggatgtgtcgggagagctgtgcactgtctcgaactcattctgcagcagccggttaaaattgttaataagcagagattataaatcattctgcagcttccctgcactgaattctcttga
ttgcagaagaaagtgtaagctctcaagtggcgttttgaacctgctgaatgttctgaaagaaactctgcgcgggaattatgcctgtgcacagagctctcttcattcatgctccctttacatgaatgaagaaacagagatctctcgaacccctgccccataa

More info: Please contact us for any question: contact@proteogenix.fr