



Reference: PTXBC007682

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MICCAP; STAM-binding protein; associated molecule with the SH3 domain of STAM; endosome-associated ubiquitin isopeptidase; testicular secretory protein Li 54; STAM binding protein		

For research use only.

GENE DESCRIPTION	STAM binding protein
GENE ID	10617
ORIGIN SPECIES	Human
ACCESSIONS	BC007682
VECTOR	pDONR223
NCBI SYMBOL	STAMBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]