

SYNJ2BP-synaptojanin 2 binding protein Gene

Reference: PTXBC007704

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARIP2; OMP25; synaptojanin-2-binding protein; activin receptor interacting protein 5; mitochondrial outer membrane protein 25; synaptojanin 2 binding protein		

For research use only.

GENE DESCRIPTION	synaptojanin 2 binding protein
GENE ID	55333
ORIGIN SPECIES	Human
ACCESSIONS	BC007704
VECTOR	pDONR223
NCBI SYMBOL	SYNJ2BP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggagagtggtattttgtcactgaggaagagatcaatctaccagaggccctcagggtgggttcaacatcgtcgtgggacagatcagcagtgatgtctcaacgacagtggtcatctacgtcagccgcatcaagaaatggggctgcggccctggatgggcggctcca
ggagggtgataagatcctttcgtaaatggccaagacctaaagaacctgctgcaccagatgctgtagacctcttctgtaatgcaggctatgctgtgtctctgagagtgacagacaggttacaggtgcagaatggacctataggacatcagagtgaaaggggacccaagtggtattccatatt
tatgggtcgtggtgccagtgtttgccctaccatggttagcagcctgggtttcatgagataccggcaacaactttga