

SPINT2-serine peptidase inhibitor, Kunitz type, 2 Gene

Reference: PTXBC007705

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DIAR3; HAI-2; HAI2; Kop; kunitz-type protease inhibitor 2; hepatocyte growth factor activator inhibitor type 2; placental bikunin; serine protease inhibitor, Kunitz type, 2; testicular tissue protein Li 183; serine peptidase inhibitor, Kunitz type 2		

For research use only.

GENE DESCRIPTION	serine peptidase inhibitor, Kunitz type, 2
GENE ID	10653
ORIGIN SPECIES	Human
ACCESSIONS	BC007705
VECTOR	pDONR223
NCBI SYMBOL	SPINT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcgcagctgtcggcctgagcggagccggcgtttctgccctgctgggatcgtctcctctctgggctcctgcccggcggaccgagaaacgcagatccacgaactctgctgggtgtcgaaggtgtgggcagatccgggctccatgcttgggtgtacaatgtcactga
cggatcctgccagctgtttgtgtatggggctgtgacggaacagcaataaactctgaccaaggaggagtgccctcaagaaatgtgccactgtcagagaaatgccacgggtgacctggccaccagcaggaaatgcagcggatctctgtcccaagtgtcccagaagcaggattctga
agaccactccagcgatagttcaactatgaagaatactgcaccgcaacgcagtcactggccttgcctgcatcctcccacgctggtactttgacgtggagaggaactcctgcaataactcactatggaggctcggggcaataagaacagctaccgctctgaggagcctgcatgc
tcgctgttccgccagcaggagaatctcccctgcccttggctcaaggtgtgcttctggcgggctgttcgtgatgtgtgatccttctcctgggacctccatggtctacctatccgggtggcacggaggaaccaggagcgtgccctgtgcaccgtctggagctccggagatga
caaggagcagctggtgaagaacacatatgtcctgtga
```