

FAM105B-family with sequence similarity 105, member B Gene

Reference: PTXBC007706

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	family with sequence similarity 105, member B
GENE ID	90268
ORIGIN SPECIES	Human
ACCESSIONS	BC007706
VECTOR	pDONR223
NCBI SYMBOL	FAM105B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccaggctgtgggctgccccctggctgcaggaccggagctcatgctgtaccagaaaaactcataagcaatacaactggatcaagcaatgaaactggactgaaattgatgggaagaatgaggacctggtgataaaaatfaagagtccttactctgctgaggaagaagt
ggcaggcttggctgaaatgagaactgctgaagcaagacagatagctgtgatgaactatcacaaatgagcggaggaatatagcctctatgaagctgtaaaattctaagctaaacagagccattgaactataatgataaagaaaaaggaaggaagtaccattttctctgctctgt
ttgctcgggacacatcaaatgaccaggacagcttctgaggaaccctcaaccaggtgggacacactggtgcttgaacaggttgaatgtcctcttgcctatgctgtgcgccacacatccagggtaccggctccaagtacaacacggagaatcatcacagctacccccaccg
accaccaaggactggccagtggttaacgctcattgctgagacgatcggcactataacatccccgtcagagtggtgaggagaccagtctatga