

TMEM128-transmembrane protein 128 Gene

Reference: PTXBC007729

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 128		
For research use only.			

GENE DESCRIPTION	transmembrane protein 128
GENE ID	85013
ORIGIN SPECIES	Human
ACCESSIONS	BC007729
VECTOR	pDONR223
NCBI SYMBOL	TMEM128
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgccatttgggtacaaaacagaaacctccacagctgttgagaaaaaggagaaacctctccaagacttaatacattctggatttggcatccattgtgtgacctattatgttgactctttaaacccttaagaaaactccacactagcagctggtttctctgtggcagtcctt
gttgcctgtcagtttatcaattgcatttactgcatagctacctggaaatggtattgtggaattggagaatatgatgcaagtaccagccttgataccattaccactgcctcttattgcagcaggaattgtctcaacattgctttatggcatgtgtgtcgttttcaactccattgtgtgtttaccagt
ttatgggggtgtcatgtttatcacactccttgatga