

## TRIM55-tripartite motif-containing 55 Gene

Reference: PTXBC007750

|                |                                                                                                                                                                                            |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                        | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | MURF-2; RNF29; muRF2; tripartite motif-containing protein 55; muscle specific ring finger 2; muscle-specific RING finger protein 2; ring finger protein 29; tripartite motif containing 55 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tripartite motif-containing 55                                                                   |
| <b>GENE ID</b>          | 84675                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC007750                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TRIM55                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgagcgcctctgaattacaattcttttccaaagcagcagaccatggataacttagagaagcaactcatctgtccatctgcttagagatgttcacgaaacctgtgtgattctccctgtcagcacaacctgtgtgaaatgtgccagtgatatatttcaggccttaacccgtatttgc  
cacaagaggaggtaccaccatggcatcagggggccgaltccgtgcccacatctgtagacatgaagtgttttggatagacatgggtatatgtgacttcagaggaaacctgtgtgtgaaaatcatgtgacatctacaagcaggagtcaccaggccagaaaagaatccgaccagcccat  
gtgcgaggaaatgaagaggagcgcacatcaactctactgtctgaactgcgaagtacccacctgctctctgtgcaaggtgtttgtgtgcacacaaagacgtgccaggtgtgtccctcactcatgtgttccagagacagaagtctgagctcagtgatgcatcgccatcctctgtggcagcaac  
gatcgagtcagggaagtgcacagccagctggaagacacctgcaaaactatcgaggaatgttgcagaaaacagaacaagagctttgtgagaagtttgattacgtgtatggcattttggaggagaggagaatgaaatgacccaagtcatcccgaaacccaagaggagaactggaacat  
gtccgtctctgatcaaaaagtattctgatcatttggagaacgtctcaaaagtgttgaagtcagggaatcagtttatgatgagccagaatggcagtggttctgcagaatgccaaaacctgtctaaaaaatctcgaagcatcaaaaggcatttcagatggagaaaatagaacatggctatga  
gaacatgaaccacttcacagtcacactcaatagagaagaaaagataatcgtgaaattgacttttacagagaagatgaagatgaagaagaagaagggcggagaaaggagaaaaagaggagaagtgaggagagagcagtagaagtggaagaggtagaaaatgttcaaa  
cagagtttccaggagaagatgaaacccagaaaaagcttcagagctctctcaggtgagctgcaggtgtccctgtgggcactccagtttctctccagagccacctccagccctgccacctgtctgaggatgccctgtgacacagatttgatttgaggctcctccctccaggacagg  
ctgcagctccagcagtgagtgagctgattctgagccagctgccatactctctcttctgttgaactccctaaatgaatga