

NAB2-NGFI-A binding protein 2 (EGR1 binding protein 2) Gene

Reference: PTXBC007756

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MADER; NGFI-A-binding protein 2; EGR1 binding protein 2; melanoma-associated delayed early response protein; NGFI-A binding protein 2		

For research use only.

GENE DESCRIPTION	NGFI-A binding protein 2 (EGR1 binding protein 2)
GENE ID	4665
ORIGIN SPECIES	Human
ACCESSIONS	BC007756
VECTOR	pDONR223
NCBI SYMBOL	NAB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacagagcgcttccccacagccgagcagccgcccggcgaggggacagcgcccgccgacccctgcagccagactcaagccagtgcccagccatggcactgcctcgagcgtgggggagctgcagctgtaccgggtcctgcagcgcgccaacctcttctactatg
agaccttcacagcaggagggggacgacgtgcagcagctgtgtgagcgggtgaggaggagttctggagatcatggcactgtgggcatggccaccaagccccctccatgcccgcctgcagaaaggcactgagagagtgggcccaatccaggctcttcagtcacacagtg
ctgctgtcccgctccagcatcccgctctcaagatctctgagactgcgggtacccggaaggagcatgagcaatgggcatggcagccaggggaaaaaggcagcagtgcccgcagtgttagcccccaagagcccccttgaactggagagaagctatcaccactgcctggggac
ggggcaggggacccccgagctgtggcaggccggagcactccagatcggacgttggggcaggaggagaaggaggctggcagcaggtgggactgggggtgtccagaccgactggagccagagatgttacgacgctgtctgggagagctctggatggacattgcaggct
gtggggcatgtccaaggctga