



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SLC12A7-solute carrier family 12 (potassium/chloride transporters), member 7 Gene

Reference: PTXBC007760

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KCC4; solute carrier family 12 member 7; K-Cl cotransporter 4; electroneutral potassium-chloride cotransporter 4; potassium/chloride transporter KCC4; solute carrier family 12 (potassium/chloride transporter), member 7; solute carrier family 12 (potassium		

For research use only.

GENE DESCRIPTION	solute carrier family 12 (potassium/chloride transporters), member 7
GENE ID	10723
ORIGIN SPECIES	Human
ACCESSIONS	BC007760
VECTOR	pDONR223
NCBI SYMBOL	SLC12A7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr