



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMC6-transmembrane channel-like 6 Gene

Reference: PTXBC007766

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

EV1; EVER1; EVIN1; LAK-4P; transmembrane channel-like protein 6;
epidermodysplasia verruciformis 1; epidermodysplasia verruciformis protein 1;
expressed in activated T/LAK lymphocytes; protein LAK-4; transmembrane channel
like 6

For research use only.

GENE DESCRIPTION	transmembrane channel-like 6
GENE ID	11322
ORIGIN SPECIES	Human
ACCESSIONS	BC007766
VECTOR	pDONR223
NCBI SYMBOL	TMC6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctacacctttcggggagagctacacgggtggggcagcacctctgggcafcacagccatcaccgtctctgctcgtcggactacaaggctagcagcagaagcggggctcccgccctcagcaggacaaattcgcaccggcgtagaaggagctgctgccggagtggcagctgcggccacagcc
ccaggagcgtgtgcccggaggctgcggcagcggcgtgtgctggggctgtgtggctgtgtgtcgtgggacgcgcgtgggctgcgcgtggccgtccacgtctctctggaggttcattgatcagatgccagaggctgctggccaggaggagctgtgtcgtggtctgccccctgggtgtggc
ctctcaacattggggggccccctcattgtgccgtgtcctggcgccctggagccgcatgactccccgggtactggagtgctatgtggccatctgcaggagaacactatcctcaagctgcccactctctgggtgctgacactcattggctggggcgaggggcgctcgtcagccagtg
ctggagggaatttgtggccacgaactctaccctctcgtatgactctcctctacacacatttggcgaactctgtgtgaggaac