

ZNF346-zinc finger protein 346 Gene

Reference: PTXBC007775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Zfp346; zinc finger protein 346; double-stranded RNA-binding zinc finger protein JAZ; just another zinc finger protein		

For research use only.

GENE DESCRIPTION	zinc finger protein 346
GENE ID	23567
ORIGIN SPECIES	Human
ACCESSIONS	BC007775
VECTOR	pDONR223
NCBI SYMBOL	ZNF346
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtatcccgccgcccacgggtgcaggccgaggcgggagcggcgggcttacagcagctcggagttgctggaggccaggagccggggtgcgcttgaccgcgagagggcgccgcctgtgggaagccgtgtccggtgccagccggtggtagagag
gaagtggagcagatgatccagaagaacaaatgtcttcaccaacacccagtgtaagggttgcgccttgcttattctgagtcgccagaagctggcacattaccagagcaaaaaacatgccaaagtgaaagatacctagcaatccatggaatggagacattaagggggaacga
agaagctagactcagatcagaagagcagcagaagcaagacaagaaccagtgctgccccatctgaacatgaccttttctccctgtctgtggccagtcgcactacctggggaagaccacgcaaaagaacttaaaagctgaagcagcagtccttaagggtggaagccttgaccagaat
agagagatgatagaccagacaagttctgcagcctctgcatgcaactttcaacgacccctgtcatggctcaacaacattatgtgggcaagaacacagaaacaggaaccaagctcaaaactaatggcagcctatggcggtggtggacccctgtgtcactgactttcagctggaag
ggctacccctgcaaaacatgaagatagtgctgaactccatagaacagtacaaagctcatgtcagcggcttcaaacacagaaccagtcacaaaaacagtgatcatccctggccagattccaatgcaaggcaacccattcagaaagactcaaccaccttgaagactag