

CCL3L3-chemokine (C-C motif) ligand 3-like 3 Gene

Reference: PTXBC007783

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D17S1718; G0S19-2; LD78; LD78BETA; SCYA3L; SCYA3L1; C-C motif chemokine 3-like 1; G0/G1 switch regulatory protein 19-2; LD78-beta(1-70); chemokine (C-C motif) ligand 3-like 3; chemokine (C-C motif) ligand 3-like, centromeric; small inducible cytokine A3-l		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 3-like 3
GENE ID	414062
ORIGIN SPECIES	Human
ACCESSIONS	BC007783
VECTOR	pDONR223
NCBI SYMBOL	CCL3L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagggtccactgctgccctggcgtcctctctgcaacalggctctctgaaccaggtcctctctgcaacacttgctgctgacacgccgaccgctgctgcttcagctacacctcccgcagattccacagaatttcagctgactcttgagacgagcagccagtgctccaagccca
gtgtcatcttctaaccaagagaggccggcaggtctgtgctgacccagtgaggagtggtggtccagaatactgctgactgacctggagccgagtgctga