

GEMIN7-gem (nuclear organelle) associated protein 7 Gene

Reference: PTXBC007793

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIP3; gem-associated protein 7; gemin-7; gem nuclear organelle associated protein 7		
For research use only.			

GENE DESCRIPTION	gem (nuclear organelle) associated protein 7
GENE ID	79760
ORIGIN SPECIES	Human
ACCESSIONS	BC007793
VECTOR	pDONR223
NCBI SYMBOL	GEMIN7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaacccagtgaaacattccctgctgtgctccggtgccccggggccctgatggcttcagccgtggctttgccctgatggacgcagagcccccttgaggccagaggttctgaaatccaggagtgctccatagctcaagaatccctggaatccaggagcagcgggcacgagc
cgcccttcgggagcgttacctccgcagcctgctggccatgggtgggtcatcagggtgagcttcacgttcacgaggggtgtgctgtgtggccgccactttggagccaccgacctggatgtggccaacttctacgtgtcacagctgcagactcccatagggtgtgcaagcagaggcgtgctccg
atgtatgacattatttcatactacattcaagccataa