



Reference: PTXBC007803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HAYOS; ATPase family AAA domain-containing protein 3A; ATPase family, AAA domain containing 3A		

For research use only.

GENE DESCRIPTION	ATPase family, AAA domain containing 3A
GENE ID	55210
ORIGIN SPECIES	Human
ACCESSIONS	BC007803
VECTOR	pDONR223
NCBI SYMBOL	ATAD3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]