

FLJ22795-hypothetical protein FLJ22795 Gene

Reference: PTXBC007811

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	golgin subfamily A member 2-like		
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ22795
GENE ID	80154
ORIGIN SPECIES	Human
ACCESSIONS	BC007811
VECTOR	pDONR223
NCBI SYMBOL	FLJ22795
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgaacatcccaggggacctggagagccgggagggccatgggtggcattttcaactccgctggagccagtgcccaggaggaaacaaaggggtgtgctgccagccctggctcaccagtgccctgtgccagaagaagccagaggtagcggccccagccccagagtggggggtg
agtcgtgtttggggagaccaccgggcccgtgcaggggcatggagaagctgcagcgactttatggaaggagaaggtggacctga