

## FAM83A-family with sequence similarity 83, member A Gene

Reference: PTXBC007828

|                |                                                                                                                                                  |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                              | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                             | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | protein FAM83A; BJ-TSA-9; tumor antigen BJ-TSA-9; tumor-specific gene expressed in prostate protein; family with sequence similarity 83 member A |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | family with sequence similarity 83, member A                                                      |
| <b>GENE ID</b>          | 84985                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC007828                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FAM83A                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgagccgggtcaagcacctgggcaaaatccggaagcgtctggaagatgtcaagagccagtggtccggcagcagggctgacttagtgacaacgagagtgccggctggccacggagccctcttgatgggggttctgaagcctactggcgggtgtcagcaggaaggcga  
gggtgacttctgtctcgtggtgagggccagtagatccagggccagggcagggcgtgtgtcccccagacacctgtggagggcggaagcagggccctaaaggagtgactccagctccctacagtcggcactacttccctgtggcctcagaggcagcagccggccctac  
tgcacagctgggctcagctgagaagccctacctgaaggaaaaatccagcggccactgtgtacttccagaccgtcaagcacaacaacatcagagacctcgtccgctgcatcacccgactagccagaacatttcacccggagtggtggaaggagagatactgtgccaagtcaggc  
aggaaattcactggccaaatcgggagaaagttcatctcggactggagattgtctgtgtgtgactttacagcttcacttggtctgtcggacacgtgcaccggaacatctctccaagttcacaggccagggcgggtggagctgtttgacgaggagttccgccactctacgcctctccaag  
cctgtgatgggctgaagtccccgcggctgtcgtcccccgtcccgccggagcagcccgccaatggccgcttagcagcagcagtggtcctccagtgaccgcacgtctccaaccccttcagcggcgtcgtggcagcagccacccgggtacccgaactgacggctga