



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MDFI-MyoD family inhibitor Gene

Reference: PTXBC007836

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	I-MF; I-mfa; myoD family inhibitor; inhibitor of MyoD family a; myogenic repressor I-mf		

For research use only.

GENE DESCRIPTION	MyoD family inhibitor
GENE ID	4188
ORIGIN SPECIES	Human
ACCESSIONS	BC007836
VECTOR	pDONR223
NCBI SYMBOL	MDFI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]