

ARHGAP19-Rho GTPase activating protein 19 Gene

Reference: PTXBC007848

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rho GTPase-activating protein 19; rho-type GTPase-activating protein 19; Rho GTPase activating protein 19		

For research use only.

GENE DESCRIPTION	Rho GTPase activating protein 19
GENE ID	84986
ORIGIN SPECIES	Human
ACCESSIONS	BC007848
VECTOR	pDONR223
NCBI SYMBOL	ARHGAP19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatagactttgcatgggaccaaataatctggtgtgaaaaaccctgctgtagtagtgaagagcatcaggagatactgactgtacctgagggccaaaacaggagcaggttaa