

CC2D1B-coiled-coil and C2 domain containing 1B Gene

Reference: PTXBC007912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil and C2 domain-containing protein 1B; FRE under dual repression-binding protein 2; five prime repressor element under dual repression-binding protein 2; freud-2; coiled-coil and C2 domain containing 1B		

For research use only.

GENE DESCRIPTION	coiled-coil and C2 domain containing 1B
GENE ID	200014
ORIGIN SPECIES	Human
ACCESSIONS	BC007912
VECTOR	pDONR223
NCBI SYMBOL	CC2D1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatctgatcattgtccgggaatgaacctcccagccctccaggggtgactcccgatgacctggatgcttttgcgggttgagtttctactaccctaactcggaccaggctcaaaaagcaaaacagctgtggtgaagaacacaaactctcagaatttgatcaactctcaactaaacat
caaccgaaccacggggcttcaagagggtgatccagagcaaaaggcatcaagtttgagatctccacaaagggtccttctcagaagcgacaagctggttggcagcacacctgaaactggagcggctggagaatgagtgatgagatcagagaattgtggaggtcctggatggaagg
aagcccaccgggggaagctggaggtgaaggtgaggctgcgggagcctctgagtgccaggatgtgcagatggtcactgagaactggctggtctggagcccaggggcttgtga