



Reference: PTXBC007914

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SDPIII; protein kinase C and casein kinase substrate in neurons protein 3; SH3 domain-containing protein 6511; syndapin III; protein kinase C and casein kinase substrate in neurons 3		

For research use only.

GENE DESCRIPTION	protein kinase C and casein kinase substrate in neurons 3
GENE ID	29763
ORIGIN SPECIES	Human
ACCESSIONS	BC007914
VECTOR	pDONR223
NCBI SYMBOL	PACSIN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]