

EID1-EP300 interacting inhibitor of differentiation 1 Gene

Reference: PTXBC007944

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf3; CRI1; EID-1; IRO45620; PNAS-22; PTD014; RBP21; EP300-interacting inhibitor of differentiation 1; 21 kDa pRb-associated protein; CREBBP/EP300 inhibitor 1; CREBBP/EP300 inhibitory protein 1; E1A-like inhibitor of differentiation 1; NB4 apoptosis re		

For research use only.

GENE DESCRIPTION	EP300 interacting inhibitor of differentiation 1
GENE ID	23741
ORIGIN SPECIES	Human
ACCESSIONS	BC007944
VECTOR	pDONR223
NCBI SYMBOL	EID1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggaaatggctgagttgtccgagctgtatgaagagagcagtgacctgcagatggatgtgatgcctggcaggggtgaccttccgcagatggaggttaggcagcgggagccgggagctatccctgcgtccctcccgcagcggggcccaacagctcgaggaggaaggcccaatgga
ggagagggagggcccaatggcggcgccagaggggaaacggagccttgcctaacgggcccacgctggggagcagccagccaggtggcgggcccagacttcgagagcgaggacgagggcgaggaatttgatgactgggaggacgactacgactatcccgaagaggagc
agctcagtggtgcccgtacagagtatcagccgctcttgaagaagccgacaagatgtttctgagaacaagagaaccagccctggatggcgggttcagatgcattatgagaagacccgtttgatcagtttagctttatcgaagagctttttcactgatggtgtcaatcgtctgaccgaagaa
ctcgctgtgatgagattattgatagagtag