

MGC16025-hypothetical protein MGC16025 Gene

Reference: PTXBC008026

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC85009		

For research use only.

GENE DESCRIPTION	hypothetical protein MGC16025
GENE ID	85009
ORIGIN SPECIES	Human
ACCESSIONS	BC008026
VECTOR	pDONR223
NCBI SYMBOL	MGC16025
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgtagggagacgcgcaccgagctggagacacagatctgcgaggagtaactgcctctctaccagcactgccagggcgcatcagcgcagcaccttcagcatgcgctaactgctggagaagctggagtgctggaaggagccctgcctgcctgcctcctcctcagtt
gctggccttgccagctctagggta