



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KLK5-kallikrein-related peptidase 5 Gene

Reference: PTXBC008036

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KLK-L2; KLKL2; SCTE; kallikrein-5; kallikrein-like protein 2; stratum corneum tryptic enzyme; kallikrein related peptidase 5		

For research use only.

GENE DESCRIPTION	kallikrein-related peptidase 5
GENE ID	25818
ORIGIN SPECIES	Human
ACCESSIONS	BC008036
VECTOR	pDONR223
NCBI SYMBOL	KLK5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctacagcaagacccccctggatggtggctctgtgctgcatcacagccttgctctgggggtcacagagcattgtcttccgccaaacaatgatgttctctgtgaccaccctctaacaccgtgccccttgggagcaaccaggacctgggagctggggccgggggaagaccgccggctg
gatgacagcagcagccgcatcatcaatggaatccgactgcgatatgcacaccagccgttgaggcggccgctgttctaaggcccaaccagctctactctggggcggtgttggctcatccacagtggtgctctacggccgccactgcaggaaagaattttcagagtccgtctcggcca
ctactccctgtcaccagtattatgaatctgggcagcagatgttccaggggggtcaaatcatccccaccctggtctactcccccttgcgcacttaacgacctatgctctacaactgaacagaaagtctgccactaaagatgtgcagaccatcaactctctctcattgtccctctctgctgg
gacaaaagctgtctgtctgtctgggggacaaaagaagcccccaagtgcaactctcctaagctctcctgactgtaataacagctgctaagtgcagaaaaggtgcaggagtctctaccgccagagacatagatgacacacattctcgcggctgcagaaagcagtagagactctcgca
cgtggaatctggggcgccctgtgtctgtcgaatggctccctcagaggacactctgtctcggcggaattacccttgcggccgcccaacagacagcgtgtgtctacacgaacctctcgaattaccacagtgatccaggaacacatccctgcggccactctga