



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM41A-transmembrane protein 41A Gene

Reference: PTXBC008043

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

2900010K02Rik; transmembrane protein 41A

For research use only.

GENE DESCRIPTION	transmembrane protein 41A
GENE ID	90407
ORIGIN SPECIES	Human
ACCESSIONS	BC008043
VECTOR	pDONR223
NCBI SYMBOL	TMEM41A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcgcgcttctcggcctcttctgcttctcggcgcgcaccttgccttctacttctgctgcagcgcgactgccccgcggcgaggagctggcgctccaccgaggaggctggaggcaggctgctgtgttccctccgcacctggcagagctgcgggagctctctgaggctctctgagagt
accgggaaggagcaccaggcctacgtgttctctctctctcggcgctctacctacatcaaacaggggcttggcatccccggctccagctctctgaagtgttttagctgtgtgctgtgttggggcgaaggctggggctctctctgtgtgtgtgtgtgacctcgtgggtggcgacatgctgctacctctct
cagatatttttggcaaacagttgtgtgtgtctacttctctgataaaggctggcgctgcgcagagaagaagctggaggagacaacagaaacagctgtgttttttctattgttttggagacttttcccatgacacaaacatgggtctggaacctctcggcccaattctgaacattccatcgtgcagcttcttct
ctcagcttcttaccgtgttggatcccatataatttctgtgtgcagcagggttcatctgtcaacctcaactctctgatgtctgtggagaccttcttaagctgttggccattgccaatggtgtgcattatctctgggaacctcattaaaaattatgcagaaacatctgcaattgaatgaa
acaagtactctctaatcatctatcatcacagttaacaagaacagata