

DIRAS2-DIRAS family, GTP-binding RAS-like 2 Gene

Reference: PTXBC008065

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Di-Ras2; GTP-binding protein Di-Ras2; DIRAS family, GTP-binding RAS-like 2; distinct subgroup of the Ras family member 2; DIRAS family GTPase 2		

For research use only.

GENE DESCRIPTION	DIRAS family, GTP-binding RAS-like 2
GENE ID	54769
ORIGIN SPECIES	Human
ACCESSIONS	BC008065
VECTOR	pDONR223
NCBI SYMBOL	DIRAS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgagcagagtaacgattaccgggtggcctgtttggggctggcgtgttggcaagagctccctggtgttgaggttgtgaaaggacattccgggagagctacatcccgacgttggaagacacctaccggcaagtgatcagctgtgacaagagcatatgcacattgcagatcacc
gacacgacggggagccaccagttcccggccatgcagcggctgtccatctccaaaggcagccttcctgtgtactccattaccagccgacagtccttgaggagctcaagcccatctacgaacaatctgcgagatcaaaaggagacgtggagagcatccccatcatgctgtggg
gaacaagtgtgatgagagcccgagcgggtgagagcagcagggcggagggccttggcccgacatggaagtgtgcttcctgagacctcagccaagctcaaccataacgtgaaggagctttccaggagctgctcaacctggagaagcgcaggacgtgagctccagatcg
acgggaaaaagagcagcagcaaaaaggaaagagaagctcaaggcaagtgcgtgatcatgtga