



Reference: PTXBC008207

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H/ACA ribonucleoprotein complex non-core subunit NAF1; nuclear assembly factor 1 homolog; nuclear assembly factor 1 ribonucleoprotein		

For research use only.

GENE DESCRIPTION	nuclear assembly factor 1 homolog (S. cerevisiae)
GENE ID	92345
ORIGIN SPECIES	Human
ACCESSIONS	BC008207
VECTOR	pDONR223
NCBI SYMBOL	NAF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]