

SERF2-small EDRK-rich factor 2 Gene

Reference: PTXBC008214

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4F5REL; FAM2C; H4F5REL; HsT17089; small EDRK-rich factor 2; gastric cancer-related protein VRG107; protein 4F5-related		

For research use only.

GENE DESCRIPTION	small EDRK-rich factor 2
GENE ID	10169
ORIGIN SPECIES	Human
ACCESSIONS	BC008214
VECTOR	pDONR223
NCBI SYMBOL	SERF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgggcttctgctgccgccgcaagcagagggactcggagatcatgcagcagaagcagaaaaaggcaaacgagaagaaggaggaacccaagtagctttgtggcttcgtgtccaacctcttgcccttcgctgtgtgctggagccagtcaccacagctcgctttctcctgta
gtgctcacagggtccagcaccgatggcattcccttgccctgagctcgcagcgggtccctttgtgtcttccctcaggtagctctctccccctgggccaactccgggggtga