

CCK-cholecystokinin Gene

Reference: PTXBC008283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cholecystokinin triacontatriapeptide; prepro-cholecystokinin		
For research use only.			

GENE DESCRIPTION	cholecystokinin
GENE ID	885
ORIGIN SPECIES	Human
ACCESSIONS	BC008283
VECTOR	pDONR223
NCBI SYMBOL	CCK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacagcggcgtgtgcctgtgctgtgatggcggtactggcggtgacgcagccgtgctcccgagatcccggggctccggctgagcgggcagaggaggcgcccgtaggcagctgagggatcgagagaacggatggcgagtcggagcgacac
ggcgccctgctggcaagatacatccagcaggccggaaagctctcttgacgaatgtccatcgttaagaacctgcagaacctgacccagccacaggataagtgaccgggactacatggcctggatggatttggccgtcgcagtgccgaggagtagtagtaccctcctag