

NOL6-nucleolar protein family 6 (RNA-associated) Gene

Reference: PTXBC008298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NRAP; UTP22; bA311H10.1; nucleolar protein 6; nucleolar RNA-associated protein; nucleolar protein 6 (RNA-associated); nucleolar protein family 6 (RNA-associated)		

For research use only.

GENE DESCRIPTION	nucleolar protein family 6 (RNA-associated)
GENE ID	65083
ORIGIN SPECIES	Human
ACCESSIONS	BC008298
VECTOR	pDONR223
NCBI SYMBOL	NOL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcattgttaccaccaagaccgcaaaaactctgtgtggacacagatggaccctcagccagatcctgcagcagcttgggtcctggcagctgaagccctgcccatttagagaagcagctcatggatccccggggacgtgggacatcaggacagtggtccggccgcttggac
atttacgacgtgctgattcgcctgtctcctcgcataatccggggcaccgcccaggtgtggactcgcagctgcctcctctctccggggcctgctcagccagccggggccctcatccctgatgccctgctgggctatgacccctcagctctatctgacgagctcagggaggccttggg
gatctggccctttctctatgaccagcatgggtggagaggtgattggtgtcctctggaagcccaccagcttccagccgagccctcaaggcctcagcacaagggggcgcagtggtgatgtctcaggtggggagctagtaatggtgcccaatggtgaagcaatcctggaggacttctgtg
ctgggtgaaggcctggtgcagactgtggaggccgaagtgaaggtggactgtgtga