

C14orf139-chromosome 14 open reading frame 139 Gene

Reference: PTXBC008299

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 139
GENE ID	79686
ORIGIN SPECIES	Human
ACCESSIONS	BC008299
VECTOR	pDONR223
NCBI SYMBOL	C14orf139
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttcagcaggatgttcaggccctaaagcacagtcagcaggactcagcatcctccagcactgctctcgggtgtcatggttaacatcatcccaaccaaccacctgtccagccgagagacagcaatcataaggaggacctcgggttccccgaggatcctgggtctcttctgaaa
cgcttctctctgagctcagcaaccagggaacaccaggccagccatcccccagcactctgtggagatgagggaacaaagtcctacagtcctcttctgttctgatgagaaaggagggaagaaacatacccgagcgctgcaatatggatgacacttcaaaaagcctgtgctatgga
gtcatgatcagaaaccagatgtgtggagagggtcagcagcctgccicagagcagccagctaggcggggagtggttaatttgggactgtaccagcgatgactggctccgagccagtgctccactctatggaatgttcctgggctcagttgcttcttcttgcaggccgaggctg
ctgccactctggcagctggtagttag