

COMMD8-COMM domain containing 8 Gene

Reference: PTXBC008371

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COMM domain-containing protein 8; COMM domain containing 8		

For research use only.

GENE DESCRIPTION	COMM domain containing 8
GENE ID	54951
ORIGIN SPECIES	Human
ACCESSIONS	BC008371
VECTOR	pDONR223
NCBI SYMBOL	COMMD8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagccggaagaggggacgcccttggtggcgtgcagaagctgccggcgagctggcgccgcagctcttcacaaaataattgatggcatttggtgctgagctctgtgtaccaagattatcacactgtttgggaatcagaagaatggatgcacgttttagaatattgccaatttt
caaaagccatagttggtaaaaacttacctgatgaagagatattcagcagttgaatcagttgaattcacttcatcaagaaactatcatgaaatgcgtgaaaaagtagaaagatgaaatcaaacaggctctgtcaagagaaatagttgctatttctctgcacagctacaggtttgattggcaggta
aagcttgcactttccagtgacaagattgctgcattacgaatgccacttttaagcctgcacttagatgtaaaagaaaatggtgaagttaaacccttattctattgaaatgagtagagaggctgcagaatctaatacagtccttgggaagcagcgaataaggtgctcctgcagttgaaataa