

MED21-mediator complex subunit 21 Gene

Reference: PTXBC008380

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SRB7; SURB7; hSrb7; mediator of RNA polymerase II transcription subunit 21; RNA polymerase II holoenzyme component SRB7; RNAPII complex component SRB7; SRB7 suppressor of RNA polymerase B homolog; mediator complex subunit 21		

For research use only.

GENE DESCRIPTION	mediator complex subunit 21
GENE ID	9412
ORIGIN SPECIES	Human
ACCESSIONS	BC008380
VECTOR	pDONR223
NCBI SYMBOL	MED21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggatcggtcagcagcttcaggacgtgtgaattcgttcagatcagttttgtaatgccattggagtattgcagcaatgtgtgtcctcttcaataattcagacagcaattacaagaccagcagctaacctacagaagagtagtcccagcttttgcagcactgattg
cacgaacagcaaaagacattgatgttttgatagattcttaccagtgagaatctacagctgctttacaggctgtagctgtataagctagaagaagaaacctgaagctgctacatgtctggaggtgtgtttatcaggagacatgcttctggagaagatacaagcgcacttgcgtgata
ttgcacagtcacagctgaagacaagaagtggtaacctagccagctcttcacagactcatag