

PDE4D-phosphodiesterase 4D, cAMP-specific (phosphodiesterase E3 dunce homolog, Drosophila) Gene

Reference: PTXBC008390

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACRDYS2; DPDE3; HSPDE4D; PDE43; PDE4DN2; STRK1; cAMP-specific 3',5'-cyclic phosphodiesterase 4D; cAMP-specific phosphodiesterase PDE4D6; phosphodiesterase 4D, cAMP-specific (phosphodiesterase E3 dunce homolog, Drosophila); testicular tissue protein Li 136		

For research use only.

GENE DESCRIPTION	phosphodiesterase 4D, cAMP-specific (phosphodiesterase E3 dunce homolog, Drosophila)
GENE ID	5144
ORIGIN SPECIES	Human
ACCESSIONS	BC008390
VECTOR	pDONR223
NCBI SYMBOL	PDE4D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcagcagacaagccggacactttaacagctacgtgaagtggataatccgcatgtccaaaccgtggctgaacgaagacctgtgaaatcctgcgagaaaacctgtgcagcatgagaagccaagacagcgaggaaatcggtttctccaaagctctccagtgatctctccgag
aaattccccaggtctctgcgcagaatgcttctcagcagcaacatccccaaacagcggtttcacgggtggcacatacatgtttgatgtggacaatggcacatctcgggacggagtccttgatcccatgaccagcccagatccgggtaattctccaagcaaatgttccacagtca
acgacgggagtccttctgtatcgatccgacagcgattatgacctctctccaaagtctatgtcccggaactcctcattgccagtgatatacacggagatgacttgattgtgactcattgtctcaggtcttggccagctgcgaactgtacgaacaactttgctgattaactaattgcaagatc
gagcacctagcaaaagatcacctatgtgcaaccaaccatcatcaaaaagccaccataacaggatcctggatggaattgaacctatacgttacttgacatgtga