

## PSG4-pregnancy specific beta-1-glycoprotein 4 Gene

Reference: PTXBC008405

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PSBG-4; PSBG-9; PSG9; pregnancy-specific beta-1-glycoprotein 4; PS-beta-G-4; PS-beta-G-9; Pregnancy-specific beta-1-glycoprotein-4; pregnancy-specific beta-1-glycoprotein 9; pregnancy-specific glycoprotein 4; pregnancy-specific glycoprotein 9; pregnancy s |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | pregnancy specific beta-1-glycoprotein 4                                                         |
| <b>GENE ID</b>          | 5672                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC008405                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | PSG4                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgggaacctctcagccctccctgcacacagcgcatcaaatggaaggggctctgctcagcagcatctttaaactctggaacctgcccaccactgcccagtcagattgaagccgagccaacaaagtctccgaggggaagatgtctctactgtccacaattgccccagaat  
cttaccggctacatctggtacaaaggcacaatgagggacctctaccattacattacatcatatgtagtagacgggtgaataaattatataatgggctgcatatggtggacgagaacagcattccaatgcacccctgctgatccagaatgtcaccgggagagcagggatcctacacctac  
acatcataaaggagatgatggactagaggagtaactggacgtttcaccttcaccttacacctggagactcctaagccctccatctccagcagcaacttaaatccaggaggacatggagcctgtgagcttaacctgtgacctgagactccagacgaagctacctgtgtggaatgaat  
ggtcagagcctccctatgactcagacttgaagctgtccgaaacacacaggacctcttctatttgggtgtcacaagatatactgcagggacctatgaatgtgaataacggaaccagtgagtgccagccgagtgaccagtcacctgaatctctccaaagctgtccaagccctacac  
acaatcaacaacttaaacccagagagaataagatgtcttaaccttcacctgtgaacctaaagagtgagaactacacctacatttgggtgctaaatgtcagagcctccctgtcagtcaccgggtaaagcgacccattgaaaacaggaatcctcattctaccaatgtcagagaaatgaaaca  
ggaccttatcaatgtgaaatacgggaccgatatggtgcatccgagtgaccagtcacctgaatgtctctatgttccagacctcccagcatttacccttcattcactattaccgttcaggagaaaacctctactgtctctcctcgcagagtcacccacgggcacaaatattctggaca  
attaatgggaagtttcagctatcaggacaaaagctctctatccccaaataactacaagcatagtggtgctctatgctgtctgttaactcagccactgccaaggaagctccaatccatcacagtcaggctctgactggatattaccctga