



Reference: PTXBC008412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BLOV1; solute carrier family 35 member E3; bladder cancer overexpressed protein; bladder cancer-overexpressed gene 1 protein; solute carrier family 35, member E2		

For research use only.

GENE DESCRIPTION	solute carrier family 35, member E3
GENE ID	55508
ORIGIN SPECIES	Human
ACCESSIONS	BC008412
VECTOR	pDONR223
NCBI SYMBOL	SLC35E3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]