

## TM4SF1-transmembrane 4 L six family member 1 Gene

Reference: PTXBC008442

|                |                                                                                                                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | H-L6; M3S1; TAAL6; transmembrane 4 L6 family member 1; membrane component, chromosome 3, surface marker 1; transmembrane 4 superfamily member 1; tumor-associated antigen L6; transmembrane 4 L six family member 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | transmembrane 4 L six family member 1                                                             |
| <b>GENE ID</b>          | 4071                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC008442                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TM4SF1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtgctatgggaagtgtgcacgatgcacgtgacattctctgttggggctcgcctctctgtgcacgcggctaataatttgccttacttcccaatggggaacaaagtatgcctccgaaaccacctcagccgttcgtgtgttcttcttggcatcgtaggaggtggcctgctgatgctcctgc  
cagcatttgcctcattgggctggaacaggatgactgctgtggctgctgtggccatgaaactgtggcaaacgatgtgcgatcttctctgtattggctgctctcattggaattgcaggatcgtgctactgtgctcattgtggcagcccttggccttagcagaaggaccactatcttgcctcctg  
gccagtggaaactacaccttggcagcaccgagggccagtaccttctggatacctccacatggtccgagtgactgaaccaagcacatttgggaatggaatgtatctctgttttctatcctcttggctcttgggtgaattgaattcatcttgtgcttattcaagtaataatggagtgcttggaggca  
tatgtggcttctgtctcaccacagcaatatgactgctaa