

SUMO2-SMT3 suppressor of mif two 3 homolog 2 (*S. cerevisiae*) Gene

Reference: PTXBC008450

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSMT3; SMT3B; SMT3H2; SUMO3; Smt3A; small ubiquitin-related modifier 2; SMT3 homolog 2; SMT3 suppressor of mif two 3 homolog 2; sentrin 2; ubiquitin-like protein SMT3A; ubiquitin-like protein SMT3B; small ubiquitin-like modifier 2		

For research use only.

GENE DESCRIPTION	SMT3 suppressor of mif two 3 homolog 2 (<i>S. cerevisiae</i>)
GENE ID	6613
ORIGIN SPECIES	Human
ACCESSIONS	BC008450
VECTOR	pDONR223
NCBI SYMBOL	SUMO2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgacgaaagcccaaggaggtcaagactgagaacaacaatcatattgaagggtggcgggcaggatgttctgtgtgcagtttaagattaaggagcagacacacacttagtaactaatgaaagcctattgtgaacgacagggtgtcaatgaggcagatcagattcc
gatttgacgggcaaccaatcaatgaaacagacacactgcacagttggaaatggaggatgaagatacaattgatgtttccaacagcagacgggaggtgtctactga