



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DPM1-dolichyl-phosphate mannosyltransferase polypeptide 1, catalytic subunit Gene

Reference: PTXBC008466

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDGIE; MPDS; dolichol-phosphate mannosyltransferase subunit 1; DPM synthase complex, catalytic subunit; DPM synthase subunit 1; MPD synthase subunit 1; dolichol monophosphate mannose synthase; dolichol-phosphate mannose synthase subunit 1; dolichyl-phosph		

For research use only.

GENE DESCRIPTION	dolichyl-phosphate mannosyltransferase polypeptide 1, catalytic subunit
GENE ID	8813
ORIGIN SPECIES	Human
ACCESSIONS	BC008466
VECTOR	pDONR223
NCBI SYMBOL	DPM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctccttggaggctcagctcgtgctcctgcaggctcctggcgagctgaggaagctgcagctccacagacagaaacaattatcggctctttacctctacaaacgagcgcgagaacctgccgctcatcgtfctggctcgtfctgaaagctctccgagagtggaatcaactatgaaattat
aatcatagatgatggaagcccgatgggaacaggagatgctgctgaacatgttggagaagatcatctggctcagacagaattcttaagaccacgagagaaaaagttggcactaggaaactcatalattatcgaatgaaacatgccacaggaaactacatcattattatgatgctgatctc
caccatcaaaatttatctgaattattaggaagcaaaaggaggtaatttgatattcctcggaaactcgctacaagaagaaatgagggtgatfctggctggcgaattgaaaagaaaaataacgccgtggggccaatttttaacacagatctctgacgacagcaggagcatctgattaacagg
aagtttcagattataccgaaaagagcttctacgaaatataagaaaaatgtgttctaaaggctactcttccctgagatgagatgttctggcgcaagacagatgaattatactattggcagggttccaatatcatttgggatcgtgttatgtggaatccaagttggggaggaatgaaatagatc
ttcttcaagagattatgaactcttttctactacataa

More info: Please contact us for any question: contact@proteogenix.fr