

C11orf52-chromosome 11 open reading frame 52 Gene

Reference: PTXBC008501

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 24 open reading frame, human C11orf52; C24H11orf52		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 52
GENE ID	91894
ORIGIN SPECIES	Human
ACCESSIONS	BC008501
VECTOR	pDONR223
NCBI SYMBOL	C11orf52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacacgggtctgctgcggaggaagctggagctgcccataactttccagaagaaaaagaaagaggaagccaaacaaagggacactgaagccgcagccacaacagctgcagcagaatctccaaaggccatgaaacaacaggacatactatgaacgggtgttacacagcaagggtctcaagagagagtcagccagcctcatgtcgaagacagcaacttacattatctgacattcaagtgtgcagccgtcccatgccgggaagtgaacacgtgcatttagaaaacgctacagagatgcgacccctcgtctcccccaggccacacctgctatgacagcaagaacgggacacctggtgtga