

C21orf7-chromosome 21 open reading frame 7 Gene

Reference: PTXBC008567

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf7; HC21ORF7; TAK1L; TAKL; TAKL-1; TAKL-2; TAKL-4; MAP3K7 C-terminal-like protein; TAK1-like protein 1; TAK1-like protein 2; TAK1-like protein 4; TGF-beta activated kinase; MAP3K7 C-terminal like		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 7
GENE ID	56911
ORIGIN SPECIES	Human
ACCESSIONS	BC008567
VECTOR	pDONR223
NCBI SYMBOL	C21orf7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcagcacagccagggtacctgtgacaagcctgtacgcatgcctttagcctcaatgacgcctcagatgataacccccctgaagactccattcttggctttccagaattagaccagcagctacagccccctgccgcttgcctatgactccgaggaatccatggaggtgttcaaacagc
actgccaaatagcagaagaataccatgagggtcaaaaaggaaatcacctgcttgagcaaaggagaaggagctcattgccaagttatgacagcagaaaaaggagaaggtggatgctgctgagctgtgttcgggaattcaggctctgacggaggagaatcgacgttgaggttggccca
gtctcaatgtgtggaacaactggagaacttcgaatacagtatcagaagaggcagggtctgtcctaa