



Reference: PTXBC008602

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD2BP1; CD2BP1L; CD2BP1S; H-PIP; PAPAS; PSTPIP; proline-serine-threonine phosphatase-interacting protein 1; CD2 antigen-binding protein 1; CD2 cytoplasmic tail-binding protein; CD2-binding protein 1; PEST phosphatase-interacting protein 1; truncated proli		

GENE DESCRIPTION	proline-serine-threonine phosphatase interacting protein 1
GENE ID	9051
ORIGIN SPECIES	Human
ACCESSIONS	BC008602
VECTOR	pDONR223
NCBI SYMBOL	PSTPIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]