

TMEM80-transmembrane protein 80 Gene

Reference: PTXBC008671

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 80		
For research use only.			

GENE DESCRIPTION	transmembrane protein 80
GENE ID	283232
ORIGIN SPECIES	Human
ACCESSIONS	BC008671
VECTOR	pDONR223
NCBI SYMBOL	TMEM80
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtttatctcagcgggaactactacgacctgtatttctcgcacgctcctgatgatcagctataaaagtcagggtgttcagctatctcaccgctacgtgctcctgatcttctctctgtttctgatgggattctagaagcagttcgggttacctggcaccaggggcaacctgacaga
ggctgagaggccgctggccgcccagcctggccctcagcgtggcaccgacctctctctgcccacttctctgtttggcaggccctagtgtgtggcgaggactggccctcagcggccacgctcctggcccttcacggcctggaggccgctctcaggtgtgtgcatcgcccttcacca
ggtag