

IL1B-interleukin 1, beta Gene

Reference: PTXBC008678

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-1; IL1-BETA; IL1F2; interleukin-1 beta; IL-1 beta; catabolin; preinterleukin 1 beta; pro-interleukin-1-beta; interleukin 1 beta		

For research use only.

GENE DESCRIPTION	interleukin 1, beta
GENE ID	3553
ORIGIN SPECIES	Human
ACCESSIONS	BC008678
VECTOR	pDONR223
NCBI SYMBOL	IL1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaagtacgtgagctcgcagtgaaatgatggcttattacagtggaatgaggatgactgttctttgaagctgatggccctaaacagatgaagtgtctctccagggacgtgacgtctgccctctggatggcgccatccagctacgaatctccgaccaccactacagcaagggttc
aggcagccgctcagttgtgtgcatggacaagctgaggaagatgctgttccctgcccacagacctccaggagaatgacgtgagcactcttctccctcatctttgaagaagaacctatcttccgacacatgggataacgagcttatgtgcacgatgcacgtgtacgatcactga
actgcacgtccgggactcacagcaaaaagcttggtgatgctgtgctcatatgaactgaaagctctccacctccaggacagggatggagcaacaagtgtgtctctcatgtcctttgtacaaggagaagaagttaatacaaaaatcctgtggccttgggctcaaggaaaagaatctg
tacctgtcctgcgtgtgaaagatgataagccactctacagctggagagtgtagatcccaaaattaccaaaagaagaagatggaaaagcgatttcttcaacaagatagaatacaataacaagctggaatttgagctgcccagttcccaactggtacatcagcactctcaagcagaaa
acatgccctgtctctgggagggaccaaaagcgccagggatataactgacttccatgcaattgtgtcttctctaa