

NEUROG1-neurogenin 1 Gene

Reference: PTXBC008687

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AKA; Math4C; NEUROD3; bHLHa6; ngn1; neurogenin-1; NGN-1; class A basic helix-loop-helix protein 6; neurogenic basic-helix-loop-helix protein; neurogenic differentiation 3; neurogenic differentiation factor 3; neurogenin 1		

For research use only.

GENE DESCRIPTION	neurogenin 1
GENE ID	4762
ORIGIN SPECIES	Human
ACCESSIONS	BC008687
VECTOR	pDONR223
NCBI SYMBOL	NEUROG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagcccgcttgagacctgcatctcgacctcgactgcgccagcagcagcgagtgacatccggcttctcaccgacgaggaagactgtgccagactccaacagcagcctcgcttggggcgcccgccggcagggcgcgcccaatctcccggcgctc
tgaggttcaggggacagggacgacgagcagggagggcgccggcgccggcgagcggtccgctccgagggcgtgctgcactgctgcgcagagccggcgctcaaggccaacgacgagcgcaaccgcatgcacaacttgaacggcgccctggacgactgcg
cagcgtgctgcctcgttccccgacgacaccaagctcaccaaaatcgagacgctgcgttcgcctacaactacatctgggctctggccgagacactgcgcctggcgatcaagggtgcccggagggcggtgccgggagcgctcctgcccgccgagtgctccctgctgcccgg
tccccaaagcccgccagcagcggagtcctggggctcaggtgcggcgccgctcccgctcttgacccagtagccagcgccctcgaagacttcacctaccgccccggcgacctgttttctcttccaagctgcccagaagactgtccacacacgcccgttttcttct
taccactag